



Rhode Island Hospital
A Lifespan Partner



Brown Medical School

ADRENAL GLANDS

This protocol is used for evaluation for and characterization of adrenal masses.

CONTRAST AGENT: None used routinely. If radiologist chooses to administer contrast, use Omniscan.

SCOUT BH

CORONAL T1 2DFLASH BH

AXIAL IN AND OUT OF PHASE

AXIAL HASTE NON FAT SAT

AXIAL HASTE FAT SAT*

PAGE RADIOLOGIST TO CHECK FOR ADMINISTRATION OF CONTRAST

(Contrast may be administered to differentiate T2 bright lesions)

Coverage should begin at the level of the dome of the liver and continue to the mid kidneys.

If evaluating for pheochromocytoma, continue AXIAL HASTE FAT SAT to level of bladder.