



Rhode Island Hospital
A Lifespan Partner



Brown Medical School

BLADDER MASS

This protocol is used for evaluation of bladder mass(es) and for follow-up evaluation of patients with known bladder cancer.

CONTRAST AGENT: Omniscan

AXIAL T1 PELVIS

AXIAL T2 FAT SAT PELVIS

AXIAL VIBE PRE GAD PELVIS

AXIAL T2 FAT SAT ABDOMEN

AXIAL VIBE PRE GAD ABDOMEN

AXIAL VIBE POST GAD PELVIS (30 sec, 60 sec, 120 sec)

AXIAL VIBE POST GAD ABDOMEN

CORONAL VIBE POST GAD LG FOV