



Rhode Island Hospital
A Lifespan Partner



Brown Medical School

UTERUS: ENDOMETRIAL CA STAGING

This protocol is used for staging of endometrial cancer.

CONTRAST AGENT: Omniscan.

AXIAL T1 2D FLASH

SAGITTAL T2

AXIAL (UTERINE LONG AXIS) OBLIQUE T2

CORONAL (UTERINE SHORT AXIS) OBLIQUE T2

SAGITTAL T1 2D FLASH FAT SAT

POST CONTRAST:

DYNAMIC SAGITTAL T1 2D FLASH FAT SAT X 3 WITH FIRST RUN AT 20 SEC

AFTER INJECTION

CORONAL T1 2D FLASH FS