



**Rhode Island Hospital**  
*A Lifespan Partner*



**Brown Medical School**

## **FEMALE PELVIS/ADNEXAL MASS**

This protocol is used for evaluation of the female pelvis and for characterization of adnexal masses.

**CONTRAST AGENT:** Omniscan

**CORONAL HASTE – LARGE FOV TO COVER AT LEAST ½ KIDNEYS**

**AXIAL T1 2D FLASH**

**AXIAL T1 2D FLASH FAT SAT**

**AXIAL T2**

**SAGITTAL T2**

**CORONAL T2 FAT SAT**

**OPTIONAL POST CONTRAST:**

**AXIAL T1 2D FLASH FAT SAT**

**SAGITTAL T1 2D FLASH FAT SAT**