



**Rhode Island Hospital**  
*A Lifespan Partner*



**Brown Medical School**

## **LIMITED FEMALE PELVIS**

This protocol is used for screening for cancer originating outside the pelvis.

**CONTRAST AGENT:** Omniscan

**CORONAL HASTE – LARGE FOV TO COVER AT LEAST ½ KIDNEYS**

**AXIAL T1 2D FLASH**

**AXIAL T1 2D FLASH FAT SAT**

**AXIAL T2**

**POST CONTRAST:**

**AXIAL T1 2D FLASH FAT SAT**