



Rhode Island Hospital
A Lifespan Partner



Brown Medical School

UTERUS ANOMALY

This protocol is applied for evaluation of uterine anomalies.

CONTRAST AGENT: Omniscan.

CORONAL HASTE – LARGE FOV TO INCLUDE AT LEAST ½ KIDNEYS

AXIAL T1 2D FLASH

AXIAL T2 2D FLASH FAT SAT

SAGITTAL T2

AXIAL (UTERINE LONG AXIS) OBLIQUE T2

CORONAL (UTERINE SHORT AXIS) T2 FAT SAT

OPTIONAL POST CONTRAST:

AXIAL T1 2D FLASH FAT SAT