

INTERNAL AUDITORY CANAL

This protocol is used for specific evaluation of the internal auditory canals.

1. LOCALIZER
2. SAG SCOUT
3. AXIAL T2 FATSAT
4. AXIAL FLAIR
5. AXIAL 2DFL HEMO
6. COR T1 PRE centered parallel to the IAC'S
7. AXIAL T1 PRE centered parallel to the IAC'S
8. DWI
9. GADO COR T1 history to pre cor t1
10. GADO AXIAL T1 history to pre ax t1
11. GADO AXIAL T1 whole brain

OPTIONAL

12. AX CISS 3D

