

BRAIN : MRA

This protocol is used for routine evaluation of the vascular supply of the brain.

- Choose one of the following MRA'S based on patient history
- For COW and VBI, angle slab parallel to the hippocampus. Place top of slab at mid-Genu and through cerebellar tonsils.
- For AVM, cover entire brain, no angle necessary.
- 4 MIPS: Entire MRA rotated, entire MRA tumbled, anterior circulation rotated, and posterior circulation rotated.

1. COW .8 thick
2. VBI 1.0 thick
3. AVM 1.25 thick

