



# **SAGITTAL SINUS**

This protocol is used for evaluation of the brain as well as for evaluation of sagittal sinus thrombosis.

## **SAGITTAL SINUS BRAIN**

1. localizer
2. SAG T1
3. AXIAL T2 FATSAT
4. AXIAL FLAIR
5. AXIAL 2DFL HEMO
6. COR T1 cover from just beyond the superior sagittal sinus and as far anterior as possible.
7. MRV SAG OBL 2D TOF oblique
  - a. Post process as spin and tumble series.
8. SAG 2D PC 5venc
9. COR 2D PC 5venc run 2 or 3 times to cover just beyond the superior sagittal sinus through the pons.
10. DWI