



BRAIN : MULTIPLE SCLEROSIS

This protocol is used for evaluation of known or suspected multiple sclerosis. Contrast used to detect active MS plaques.

1. LOCALIZER
2. SAG T1
3. SAG FLAIR
4. AXIAL PD
5. AXIAL T2 FATSAT
6. AXIAL FLAIR
7. AXIAL 2DFL HEMO
8. DWI

FOR GADO:

9. LOCALIZER
10. SAG T1
11. SAG FLAIR
12. AXIAL PD
13. AXIAL T2 FS
14. AXIAL 2DFL HEMO
15. AXIAL T1 PRE
16. GADO AXIAL FLAIR
17. GADO DWI
18. GADO AXIAL T1 run 2x 5 skip 0 at a delay of 5 minutes.