

## BRAIN AND ORBITS

This protocol is used for evaluation of the brain and orbits.

1. localizer
2. SAG T1
3. AXIAL T2 FATSAT
4. AXIAL FLAIR
5. AXIAL 2DFL HEMO
6. AXIAL T1 PRE centered parallel to the optic nerve
7. CORONAL STIR to cover from optic chiasm to globe. Only cover entire globe if pathology is seen extending into globe.
8. DWI
9. GADO AXIAL T1 FS ORBIT
10. GADO COR T1 FS ORBITS
11. GADO AXIAL T1 WHOLE BRAIN



