



PEDIATRIC BRAIN

This protocol is used for routine evaluation of the brain in patients <2 YEARS OLD.

- If gado, drop sag MPR unless development delay or new brain tumor.
- Use extremity coil if possible for better SNR.
- Children 2yrs and older receive routine brain protocol.
- Reconstruct sag MPR into axial plane using MPR tool at a thickness of 1.5.

1. LOCALIZER
2. SAG T1
3. AXIAL T2 FATSAT
4. AXIAL FLAIR
5. AXIAL 2DFL HEMO
6. DWI
7. SAG MPR

IF GADO ADD:

8. AXIAL T1 PRE
9. GADO AXIAL T1
10. GADO COR T1
11. GADO SAG T1