



BRAIN : ROUTINE

This protocol is used for routine evaluation of the brain.

- Patients with VAGAL NERVE STIMULATORS **must** be scanned on the vision unless using other transmit/receive coil (ex. Extremity coil) for brain imaging.
- Use Claustro protocol for patients that require shorter scan times.
- Use Turbo protocol for patients requiring even shorter scan times.
- For trauma, add Cor T1.

1. LOCALIZER
2. SAG T1
3. AXIAL T2 FATSAT
4. AXIAL FLAIR
5. AXIAL 2DFL HEMO
6. DWI

FOR GADO ADD:

1. AXIAL T1 PRE
2. GADO AXIAL T1
3. GADO COR T1
4. GADO SAG T1