

BRAIN : WAND

This protocol is used for evaluation of brain tumors for surgical planning.

- If patient has had a brain with contrast study within two weeks of the wand study, run Sag T1 and Sag T2.
- If the wand study is scanned with the brain gado study, run Sag T1 only.
- **DO NOT ANGLE** wand sequences.
- **DO NOT** place anything (pads, tape, etc.) next to patients head. There must be no compression of the patients head to assure accurate measurements in the OR.
- Send wand sequences to BRLAB and AESCUL.
- Contrast: MAGNEVIST

1. LOCALIZER
2. SAG T1 MPR
3. SAG T2 FS