

PITUITARY: Non-Dynamic

This protocol is used for non-dynamic evaluation of the pituitary gland. This is used for all pituitary cases other than those involving a microadenoma (usually a prolactinoma).

- Used for all MACROadenoma indications.
- Requires FULL dose of contrast.
- Contrast Agent: Magnevist
- Scanned immediately after injection.

1. LOCALIZER
2. SAG T1 centered over the pituitary
3. AXIAL T2 FATSAT entire brain
4. AXIAL FLAIR entire brain
5. AXIAL 2DFL HEMO entire brain
6. COR T1 PRE centered over the pituitary
7. GADO COR T1 history to cor T1 pre
8. GADO SAG T1 centered over the pituitary
9. GADO AXIAL T1 entire brain
10. DWI

