



GAMMA: AVM

This protocol is used for patients that are receiving gamma knife treatment for brain avm.

- **DO NOT ANGLE**
- Be sure to run gado AVM MRA immediately following contrast injection.

1. LOCALIZER
2. AXIAL T2 FATSAT
3. AXIAL FLAIR
4. AXIAL T1 PRE 5mm
5. AVM MRA
6. GADO AVM MRA no sat pulse
7. GADO AXIAL T1 5mm
8. GADO CORONAL 2mm
9. GADO AXIAL T1 2mm