



GAMMA FOLLOW-UP: AVM

This protocol is used for patients that have had a previous gamma knife treatment for arterial-venous malformations.

- Angle as you would for a routine brain study.
- Contrast: MAGNEVIST
- Cover entire brain, no angle, with MRA sequence.

1. LOCALIZER
2. SAG T1
3. AXIAL T2 FATSAT
4. AXIAL FLAIR
5. AXIAL 2DFL
6. AXIAL T1 PRE
7. DWI-DTI
8. AVM MRA PRE
9. GADO AVM MRA no sat pulse
10. AXIAL MPRAGE
11. CORONAL MPRAGE