

GAMMA: CAVERNOUS MALFORMATION

This protocol is used for patients that are receiving gamma knife treatment for brain cavernous malformations.

- **DO NOT ANGLE**
- Be sure to include ALL fiducials in FOV.

1. LOCALIZER
2. AXIAL T2 FATSAT
3. AXIAL FLAIR
4. AXIAL T1 5mm PRE cover entire brain
5. AXIAL 2DFL 2mm: cover entire brain
6. GADO AXIAL T1 MPR 2mm
7. GADO CORONAL T1 MPR 2mm