



## GAMMA FOLLOW-UP : CAVERNOUS MALFORMATION

This protocol is used for patients that are receiving follow-up to a gamma knife treatment for brain cavernous malformations.

- **Angle as you would** for a routine brain
- Contrast: Magnevist

1. LOCALIZER
2. AXIAL T2 FATSAT
3. AXIAL FLAIR
4. AXIAL T1 5mm PRE cover entire brain
5. AXIAL 2DFL 2mm: cover entire brain
6. GADO AXIAL T1 MPR 2mm
7. GADO CORONAL T1 MPR 2mm