



GAMMA: IAC

This protocol is used for patients that are receiving gamma knife treatment for vestibular schwannoma.

- **DO NOT ANGLE**
- Be sure to include ALL fiducials in FOV.

1. LOCALIZER
2. AXIAL T2 FATSAT
3. AXIAL FLAIR
4. Cor T1 scout
5. AXIAL T1 PRE 5mm
6. GADO AXIAL T1 5mm
7. GADO CORONAL T1 2mm centered over the IACs
8. GADO AXIAL T1 2mm centered over the IACs
9. GADO AXIAL 3DFL T1 1mm history to gado ax t1 2mm
10. GADO AXIAL CISS 3D 1mm history to gado ax t1 2mm