



GAMMA FOLLOW-UP: IAC

This protocol is used for patients that have had a previous gamma knife treatment for vestibular schwannoma.

- Angle as you would for a routine IAC study.
- Contrast: 1.0 mmol/kg magnevist

1. LOCALIZER
2. AX T2 FS
3. AXIAL FLAIR
4. COR T1 2mm PRE thru the IAC's
5. AX T1 2mm PRE thru the IAC's
6. DWI-ADC
7. GADO AX T1 2mm
8. GADO COR T1 2mm