



## GAMMA: MENINGIOMA

This protocol is used for patients that are receiving gamma knife treatment for brain meningiomas.

- **DO NOT ANGLE**
- Be sure to include ALL fiducials in FOV.

1. LOCALIZER
2. AXIAL T2 FATSAT
3. AXIAL FLAIR
4. AXIAL T1 PRE 5mm
5. GADO AXIAL T1 5mm
6. GADO CORONAL T1 2mm MPR
7. GADO AXIAL T1 2mm MPR

Check with Gamma Doctor to see if FLASH and CISSS sequences are needed.

8. GADO AXIAL 3DFL 1mm
9. GADO AXIAL CISS 1 mm