

GAMMA FOLLOW-UP: MENINGIOMA

This protocol is used for patients that have had a previous gamma knife treatment for meningiomas.

- Angle as you would for a routine brain study.
- Contrast: MAGNEVIST
- This protocol should be patient specific based on size and location of meningioma. Check with Radiologist for position of post gado sequences.

1. LOCALIZER
2. SAG T1
3. AXIAL T2 FATSAT
4. AXIAL FLAIR
5. AXIAL T1 PRE
6. DWI-DTI
7. GADO AXIAL T1 5 skip 0
8. GADO COR T1 2mm centered over mening.
9. GADO AXIAL T1 2mm centered over mening.