



GAMMA FOLLOW-UP: METS

This protocol is used for patients that have had a previous gamma knife treatment for brain mets.

- Angle as you would for a routine brain study.
- Contrast: double dose multihance.

1. LOCALIZER
2. SAG T1
3. AXIAL T2 FATSAT
4. AXIAL FLAIR
5. AXIAL 2DFL HEMO
6. AXIAL T1 PRE
7. DWI-DTE
8. GADO AX MPRAGE be sure to cover above scalp through the foramen magnum
9. GADO COR MPRAGE