

GAMMA: PIT

This protocol is used for patients that are receiving gamma knife treatment for pituitary adenomas.

- **DO NOT ANGLE**
- Be sure to include ALL fiducials in FOV.

1. LOCALIZER
2. AXIAL T2 FATSAT
3. AXIAL FLAIR
4. AXIAL T1 PRE 5mm
5. GADO AXIAL T1 5mm
6. GADO CORONAL T1 2mm centered over the PIT
7. GADO AXIAL T1 2mm centered over the PIT
8. GADO AXIAL 3DFL T1 1mm history to gado ax t1 2mm
9. GADO AXIAL CISS 3D 1mm history to gado ax t1 2mm