



Rhode Island Hospital
A Lifespan Partner



Brown Medical School

GAMMA FOLLOW-UP: MACRO PITUITARY

This protocol is used for patients that have had a previous gamma knife treatment for pituitary macroadenomas.

- Angle as you would for a routine pituitary study.
- Contrast: 1.0 mmol/kg magnevist

1. LOCALIZER
2. SCOUT SAG
3. AXIAL T2 FATSAT
4. AXIAL FLAIR
5. COR T1 2mm PRE
6. GADO COR T1 2mm
7. GADO AX T1 2mm
8. GADO AXIAL T1 full brain 5 skip 0