



LUMBAR SPINE: EPIDURAL ABSCESS

This protocol is used to evaluate for epidural abscess.

- **Cover from T8 through S2**
 1. Scouts
 2. SAG T1 TSE ___ IPAT
 3. AXIAL T1 STACK PRE
 4. GADO SAG T1 FS TSE ___ IPAT
 5. GADO AXIAL T1 STACK
 6. GADO STIR
 7. GADO SAG T2 TSE ___ IPAT
 8. GADO AXIAL T2 STACK Restore